



LAB LOW TEMPERATURE FREEZER -45° C

MODELS	VOLUME (Lt)	INNER DIMENSIONS (mm)			OUTER DIMENSIONS (mm)		
		W	D	H	W	D	H
LABLTF580	580 L	600	640	1500	740	824	1994

CAPACITY	580 L
TEMPERATURE RANGE	-45° C / -10° C
AMBIENT TEMPERATURE	MAX. +32° C
NO. OF COMPRESSOR	1
POWER	230V / 50 HZ
COMSUMPTION	LOW POWER CONSUMPTION
INSTALLATION	PLUG & PLAY
ALARMS	VISUAL / ACOUSTIC ALARM • OPEN DOOR ALARM • ADJUSTABLE HIGH / LOW TEMPERATURE ALARM • CONTACT REMOTE ALARM
INTERFACE	RS485
ISOLATION	70 MM
WEIGHT	125 KG
SHELVES	4
MAX. SHELF LOAD	30 KG
CONTROLLOR	ELECTRONIC MICROPROCESSOR
CASTORS LOCKABLE	4 CASTORS - 2 LOCKABLE
DEFROST	MANUAL
COOLING	DIRECT COOLING
DOOR KEY LOCK	YES
DOOR	MAGNETIC DOOR GASKET • SELF-CLOSING DOOR • SOLID DOOR
NO. OF DOORS	1
REFRIGERANT	MIX
CHAMBER MATERIAL	STAINLESS STEEL
EXTERIOR MATERIAL	POWDER COATED GALVANIZED STEEL

We reserve the right to change specifications without notice. Subject to confirmation, availability and errors.

B SCIENCE

Bavnevej 10 • 6580 Vamdrup • Denmark
 info@bscienceglobal.com • www.bscienceglobal.com
 +45 287 287 62

