

## LAB LOW TEMPERATURE FREEZER -45° C

MODELS	VOLUME (Lt)	INNER DIMENSIONS (mm)			OUTER DIMENSIONS (mm)		
		W	D	H	W	D	H
LABLTF580	580 L	600	640	1500	740	824	1994



CAPACITY	580 L
TEMPERATURE RANGE	-45° C / -10° C
AMBIENT TEMPERATURE	MAX. +32° C
NO. OF COMPRESSOR	1
POWER	230V / 50 HZ
ALARMS	VISUAL / ACOUSTIC ALARM ■ OPEN DOOR ALARM ■ ADJUSTABLE HIGH / LOW TEMPERATURE ALARM ■ VOLTAGE FREE CONTENT FOR REMOTE ALARM
ISOLATION	70 MM
WEIGHT	125 KG
SHELVES	4
MAX. SHELF LOAD	30 KG
CONTROLLOR	ELECTRONIC MICROPROCESSOR
CASTORS LOCKABLE	4 CASTORS - 2 LOCKABLE
DEFROST	MANUAL
COOLING	DIRECT COOLING
DOOR KEY LOCK	YES
DOOR	MAGNETIC DOOR GASKET ■ SELF-CLOSING DOOR ■ SOLID DOOR
NO. OF DOORS	1
REFRIGERANT	MIX
CHAMBER MATERIAL	STAINLESS STEEL
EXTERIOR MATERIAL	POWDER COATED GALVANIZED STEEL

We reserve the right to change specifications without notice. Subject to confirmation, availability and errors.