

LAB REFRIGERATOR +1° / +10° C - GLASS DOOR

MODELS	VOLUME (Lt)	INNER DIMENSIONS (mm)			OUTER DIMENSIONS (mm)		
		W	D	H	W	D	H
LABREF628G	628 L	603	694	1500	740	830	2037
					SHIPPING DIMENSION (mm)		
					1000	900	2250



CAPACITY	628 L
TEMPERATURE RANGE	+1° C / +10° C
AMBIENT TEMPERATURE	MAX. + 32° C
NO. OF COMPRESSOR	1
POWER	230V / 50 HZ
ALARMS	VISUAL / ACOUSTIC ALARM ▪ OPEN DOOR ALARM ▪ ADJUSTABLE HIGH / LOW TEMPERATURE ALARM ▪ VOLTAGE FREE CONTENT FOR REMOTE ALARM
ISOLATION	70 MM
WEIGHT / SHIPPING WEIGHT	125 KG / 175 KG
SHELVES	3
MAX. SHELF LOAD	30 KG
CONTROLLOR	ELECTRONIC MICROPROCESSOR
CASTORS LOCKABLE	4 CASTORS - 2 LOCKABLE
DEFROST	AUTOMATIC
COOLING	AIR FORCED COOLING
LIGHTING	YES
DOOR KEY LOCK	YES
DOOR	MAGNETIC DOOR GASKET ▪ SELF-CLOSING DOOR ▪ GLASS DOOR
NO. OF DOORS	1
REFRIGERANT	R290
CHAMBER MATERIAL	STAINLESS STEEL
EXTERIOR MATERIAL	POWDER COATED GALVANIZED STEEL

We reserve the right to change specifications without notice. Subject to confirmation, availability and errors.