

LAB COMBI REFRIGERATOR/FREEZER +1° C / +10° C -25° C / -10° C

MODEL	VOLUME (Lt)	INNER DIMENSIONS (mm)			OUTER DIMENSIONS (mm)		
		W	D	H	W	D	H
LABCOMBI650	2 X 319 L	2 X 600	2 X 760	2 X 700	720	880	2037

GENERAL DATA

CAPACITY	2 X 319 L
TEMPERATURE RANGE	+ 1° C / +10° C & -25° C / -10° C
AMBIENT TEMPERATURE	Max. +32° C
NO. OF COMPRESSOR	2
POWER	230V / 50HZ
COMSUMPTION	LOW POWER CONSUMPTION • 2,6 KWH / 24 HOURS
ISOLATION	60MM
INTERFACE	RS485
WEIGHT	175 KG
NO. OF SHELVES	4 PIECES • MAX. LOAD 30 KG PR. SHELVE
CONTROLLOR	ELECTRONIC MICROPROCESSOR
DEFROST	AUTOMATIC
COOLING	AIR FORCED COOLING
DOOR KEY LOCK	YES
DOOR	MAGNETIC DOOR GASKET • SELF-CLOSING DOOR • SOLID DOOR
NO. OF DOORS	2
REFRIGERANT	R290
CHAMBER MATERIAL	STAINLESS STEEL
EXTERIOR MATERIAL	POWDER COATED GALVANIZED STEEL



ALARMS

- VISUAL / ACOUSTICS ALARM
- OPEN DOOR ALARM
- ADJUSTABLE HIGH / LOW TEMPERATURE ALARM
- CONTACT REMOTE ALARM