

LAB REFRIGERATOR +1° / +10° C - SOLID DOOR

MODELS	VOLUME (Lt)	INNER DIMENSIONS (mm)			OUTER DIMENSIONS (mm)		
		W	D	H	W	D	H
LABREF100	100 L	510/510	396/174	394/236	650	620	860
					SHIPPING DIMENSION (mm)		
					750	720	1000



CAPACITY	100 L
TEMPERATURE RANGE	+1° C / +10° C
AMBIENT TEMPERATURE	MAX. + 40° C
NO. OF COMPRESSOR	1
POWER	230V / 50 HZ
ALARMS	VISUAL / ACOUSTIC ALARM ▪ OPEN DOOR ALARM ▪ ADJUSTABLE HIGH / LOW TEMPERATURE ALARM ▪ VOLTAGE FREE CONTENT FOR REMOTE ALARM
ISOLATION	70 MM
WEIGHT / SHIPPING WEIGHT	145 KG / 235 KG
SHELVES	1
MAX. SHELF LOAD	30 KG
CONTROLLOR	ELECTRONIC MICROPROCESSOR
CASTORS LOCKABLE	4 FEET
DEFROST	AUTOMATIC
COOLING	DYNAMIC COOLING
DOOR KEY LOCK	YES
DOOR	MAGNETIC DOOR GASKET ▪ SELF-CLOSING DOOR ▪ SOLID DOOR
NO. OF DOORS	1
REFRIGERANT	R290
CHAMBER MATERIAL	STAINLESS STEEL
EXTERIOR MATERIAL	POWDER COATED GALVANIZED STEEL

We reserve the right to change specifications without notice. Subject to confirmation, availability and errors.