

CONSTRUCTION

✓ Internal Material Stainless Steel	✓ Replaceable Door Gasket	✓ Key Lock
✓ Self closing doors	✓ Castors	✓ Internal lighth
✓ Exterior Material Powder coated Steel	✓ Isolastion 70mm	✓ 4 Shelf

ALARMS

✓ Visual/Acoustic Alarm	✓ Open Door Alarm	✓ Contact Remote Alarm
✓ Power Failure Alarm	✓ Probe Failure Alarm	✓ Compressor Alarm
✓ Adjustable High / Low Temperature Alarm	✓ Track of Alarms & Events to be Noted Automatically	✓ Approx 24 HRS. of Battery Backup for Alarms and Logging

DISPLAY

✓ LCD Display	✓ Integrated Datalogger	✓ USB - Download
✓ Memory For Over 15 Years Of Data		

DOCUMENTATION

✓ IQ	✓ OQ
------	------

We reserve the right to change specifications without notice. Subject to confirmation, availability and errors.