

## BASIC LAB REFRIGERATOR +1/+10° C

| MODELS        | VOLUME (Lt) | INNER DIMENSIONS (mm) |     |      | OUTER DIMENSIONS (mm) |     |      | SHIPPING DIMENSION (mm) |      |      |
|---------------|-------------|-----------------------|-----|------|-----------------------|-----|------|-------------------------|------|------|
|               |             | W                     | D   | H    | W                     | D   | H    | W                       | D    | H    |
| <b>LBR700</b> | 618 L       | 540                   | 740 | 1460 | 710                   | 880 | 1998 | 800                     | 1000 | 2250 |



|                                 |                               |
|---------------------------------|-------------------------------|
| <b>CAPACITY</b>                 | 618 L                         |
| <b>TEMPERATURE RANGE</b>        | +1° C / +10° C                |
| <b>AMBIENT TEMPERATURE</b>      | MAX.+ 32° C                   |
| <b>POWER</b>                    | 230 V / 50 HZ /               |
| <b>ENERGY CONSUMPTION</b>       | 4,0 Kw / 24 Hrs.              |
| <b>WEIGHT / SHIPPING WEIGHT</b> | 160 Kg / 245 Kg               |
| <b>MAX. SHELF LOAD</b>          | 30 Kg Per. Shelf              |
| <b>COOLING TECHNOLOGY</b>       | Air Forced Cooling System     |
| <b>INTERFACE</b>                | Potential Free Remote Contact |
| <b>DOCUMENTATION</b>            | IQ • OQ                       |

### ENVIRONMENTAL EFFECTS

GWP 3 (R290)

ODP 0 (R290)



We reserve the right to change specifications without notice. Subject to confirmation, availability and errors.

# TECHNICAL DATASHEET

## CONSTRUCTION

|                                         |                           |                   |
|-----------------------------------------|---------------------------|-------------------|
| ✓ Internal Material Stainless Steel     | ✓ Replaceable Door Gasket | ✓ Key Lock        |
| ✓ Self closing doors                    | ✓ Castors                 | ✓ Internal lighth |
| ✓ Exterior Material Powder coated Steel | ✓ Isolastion 70mm         | ✓ 4 Shelf         |
| ✓ Auto Defrost                          |                           |                   |

## ALARMS

|                                           |                                                      |                                                           |
|-------------------------------------------|------------------------------------------------------|-----------------------------------------------------------|
| ✓ Visual/Acoustic Alarm                   | ✓ Open Door Alarm                                    | ✓ Contact Remote Alarm                                    |
| ✓ Power Failure Alarm                     | ✓ Probe Failure Alarm                                |                                                           |
| ✓ Adjustable High / Low Temperature Alarm | ✓ Track of Alarms & Events to be Noted Automatically | ✓ Approx 24 HRS. of Battery Backup for Alarms and Logging |

## DISPLAY

|                                    |                         |                  |
|------------------------------------|-------------------------|------------------|
| ✓ LCD Display                      | ✓ Integrated Datalogger | ✓ USB – Download |
| ✓ Memory For Over 15 Years Of Data |                         |                  |

## DOCUMENTATION

|      |      |
|------|------|
| ✓ IQ | ✓ OQ |
|------|------|

We reserve the right to change specifications without notice. Subject to confirmation, availability and errors.